



Variety Disposable gloves for your choice

一次性 PVC 手套
Disposable vinyl gloves

一次性丁腈手套
Disposable nitrile gloves





Serious work for each gloves





Specification for Vitrile Gloves



规格 SIZE	重量 Weight(g)	长度 Length(mm)	手掌宽度 Palm width(mm)	厚度 Thickness(mm)		
				指尖 Finger lip	掌部 Palm	腕部 Cuff
XS	4.0±0.2	235 ±5	75±5	0.07±0.02	0.06±0.02	0.06 ±0.02
S	4.5±0.2	235 ±5	85±5	0.07±0.02	0.06±0.02	0.06±0.02
M	5.0±0.2	235 ±5	95±5	0.07 ±0.02	0.06±0.02	0.06 ±0,02
L	5.5±0.2	235 ±5	105±5	0.07±0.02	0.06±0.02	0.06±0.02
XL	6.0±0.2	235 ±5	115±5	0.07±0.02	0.06±0.02	0.06 ±0.02
拉伸强度 TENSILE STRENGTH				≥ 13Mpa		
伸长率 ELONGATION				» 330%		
证书 CERTIFICATES				CE, FDA, ISO9001, ISO13485, SGS		
质量标准 QUALITY STANDARDS				EN455,EN374		



Specification for Nitrile Gloves



规格 SIZE	重量 Weight(g)	长度 Length(mm)	手掌宽度 Palm width(mm)	厚度 Thickness(mm)		
				指尖 Finger lip	掌部 Palm	腕部 Cuff
XS	3.0±0.2	235±5	75±5	0.07 ±0.02	0.07 ±0.02	0.06 ±0.02
S	3.5±0.2	235±5	85±5	0.07 ±0.02	0.07 ±0.02	0.06 ±0.02
M	4.0±0.2	235±5	95±5	0.07 ±0.02	0.07 ±0.02	0.06 ±0.02
L	4.5±0.2	235±5	105±5	0.07 ±0.02	0.07 ±0.02	0.06±0.02
XL	5.0±0.2	235±5	115±5	0.07 ±0.02	0.07±0.02	0.06±0.02
拉伸强度 TENSILE STRENGTH				, 20 Mpa		
伸长率 ELONGATION				» 500%		
证书 CERTIFICATES				CE, SGS		

质量标准 QUALITY STANDARDS

EN455, EN374



Nitrile Gloves

丁月青手套

产品名称 Product name: Nitrile gloves

级别 Level: AQL1.5, AQL2.5,

AQL4.0 材质 Material: Nitrile

型号 Size: S/M/L/XL

类型 Type : Powder Free/ Powdered 长

度 Length: 9 inches

表面 Surface: Fully Textured

颜色 Color: Blue



Application:

Widely used in examination, hospital, clinic, industries and other areas.



